

Beating the limitations of surface-bound scanning probe microscopy

Exploring nanoscale 3D physical properties of advanced materials and devices

Oleg Kolosov, Lancaster University, UK

Lancaster University Oleg Kolosov's group Bruker European material AFM user meeting 2021


AFM and SPM outstanding images of ... surfaces

Nanomechanics 	Surface potential (KPFM) 	Electrochemistry
Polymers 	Semiconductors 	Metals

<https://www.bruker.com/en/products-and-solutions/microscopes/materials-afm/afm-modes.html>

Lancaster University Oleg Kolosov's group Bruker European material AFM user meeting 2021

Subsurface measurements via ultrasonic force microscope (UFM) - graphite nanoflakes




Lancaster University

Oleg Kolosov's group Bruker European material AFM user meeting 2021




Looking into the 3D structures – mechanical cleavage and FIB



- Wang, A. and T. Chien (2018) Physics Letters A 382(11): 739-748.

The cleavage is a real art and suitable for epitaxies only



SEM cross-section of a GaN HEMT device.
Microelectronics Reliability 54, 1785-1789, (2014)


FIB area is small, not easily used in SPM, and Ga implantation change the properties.

Lancaster University

Oleg Kolosov's group Bruker European material AFM user meeting 2021



“Slicing” the sample directly in AFM - nanotomography and “Scalpel AFM”



Luria, J. et al, (2016) Nature Energy 1: 16150.

Chen, S. at al, Advanced Functional Materials 28(52): 1802266.


Great tool, but area is rather small, and section of harder materials like Si, iii-v is very difficult.



Oleg Kolosov's group Bruker European material AFM user meeting 2021



Can we use Ar Ion Cross-sectioning?




Low damage to the surface, but near surface layers are damaged during the preparation.



Oleg Kolosov's group Bruker European material AFM user meeting 2021



New - SPM friendly nano-cross-sectioning Beam Exit Cross-sectional Polishing (BEXP)



Method and apparatus for ion beam polishing. USA. 9,082,587.

Kolosov, O. V. et al, *Nanotechnology* 22(18): 185702.



Oleg Kolosov's group Bruker European material AFM user meeting 2021



Benefits of using BEXP nano-cross-sectioning

Thickness x10


1. Near-atomic surface roughness.
2. Low implantation rate of the Ar-ions.
3. Projection of the layer's thickness increases around 10 times.
4. Small angle from the surface (~10°), easy access with the SPM probes.



Oleg Kolosov's group Bruker European material AFM user meeting 2021



Result – BEXP based cross-sectional SPM (xSPM)



SPM-compatible near-atomically-flat section of sample subsurface layers. Thermal SiO₂ on Si with nanomechanical ultrasonic force microscopy map.


US/EU patents- US 9082587, EP 2537017 B1.



Oleg Kolosov's group Bruker European material AFM user meeting 2021



3D cross-sectioning for SPM – rivaling TEM resolution



Robson, A. J. et al *Acs Applied Materials & Interfaces* **2013**, 5 (8), 3241-3245.

xSPM vs cross-sectional TEM


Multimode SPM in both Tapping Mode AFM and nanomechanical ultrasonic force microscopy (UFM). Oxidised Al_xGa_{1-x}As layers protrude above GaAs layers by an amount which varies with Al content x, allowing identification during imaging.



Oleg Kolosov's group Bruker European material AFM user meeting 2021



xSPM of CVD layered vdW materials





CVD grown transition metal dichalcogenide (TMD) WS₂ on Si substrate.
SPM nanomechanical contrast, low defects top layer and the high defect density bottom thickness measured with nm precision.



Oleg Kolosov's group Bruker European material AFM user meeting 2021



Gate oxides



Topography

UFM




Oleg Kolosov's group Bruker European material AFM user meeting 2021



Energy storage materials - Na ion batteries

+




• NTO




Oleg Kolosov's group Bruker European material AFM user meeting 2021




KPFM (surface potential) xSPM mode - thin film solar cells (CdSe/CdS) grain structure



Topography



KPFM amplitude




KPFM phase



Oleg Kolosov's group Bruker European material AFM user meeting 2021



SThM (thermal conductivity) xSPM - Si/GeSi_x/GeSn_y multilayer MBE structure



Thermal resistance SThM map and profiles.


Inset: 3D topography overlaid with SThM response.




Oleg Kolosov's group Bruker European material AFM user meeting 2021




Devices – vertical cavity surface emitting laser (VCSEL) in xSPM - DBR and MQW structures.



VCSEL section




UFM - nanomechanics



Oleg Kolosov's group Bruker European material AFM user meeting 2021




VCSEL active area zoom in – KPFM and SSRM (conductivity) xSPM modes



Oleg Kolosov's group Bruker European material AFM user meeting 2021




Multiple quantum wells iii-v on Si structures KPFM xSPM cross-section



Oleg Kolosov's group Bruker European material AFM user meeting 2021




iii-v on Si structures multiple GaAs/AlGaAs/InGaAs quantum wells - KPFM xSPM



Oleg Kolosov's group Bruker European material AFM user meeting 2021



3D mapping of nanowires - nanoporous supercapacitor layers in Si




The structure of deep-etch vertical nanopores in the Si is observed from through the thickness of material. The structure of the nanopore-substrate interface and oxide pore clogging are clearly observed. (Sample courtesy M. Prunilla, VTT, Finland.)



Oleg Kolosov's group Bruker European material AFM user meeting 2021



Carbon nanotubes CNT thermal interface materials (SThM mode of xSPM)




CNT “forest” is “trimmed” creating gradient of CNT height via BEXP™, enabling absolute measurements of intrinsic thermal conductance of the TIM (Sample courtesy O. Bezencenet, Thales, France).



Oleg Kolosov's group Bruker European material AFM user meeting 2021



xSPM of non-planar samples - GaN NWs - embedding in Spin-On-Glass




- Deposit the Spin-On-Glass (SOG).
- Spin and bake the sample.
- BEXP sectioning.



Oleg Kolosov's group Bruker European material AFM user meeting 2021




UFM nanomechanical mapping of GaN NWs



Oleg Kolosov's group Bruker European material AFM user meeting 2021




Electron affinity (KPFM) maps of GaN NWs



Oleg Kolosov's group Bruker European material AFM user meeting 2021



Piezoelectric response GaN NWs = NW polarity



The domains of opposite polarity exist in GaN NW's. BEXP-SPM can provide insight on what trigger particular polarity growth.



Oleg Kolosov's group Bruker European material AFM user meeting 2021



EPSRC
Engineering and Physical Sciences Research Council

GRAPHENE NOWAND

RESEARCH COUNCILS UK



Consiglio Nazionale delle Ricerche

Lancaster SPM group: M. Mucientes, Y. Chen, A. Niblett, L. Forcieri, J. Speice, C. Evangelisti, E. Castanon,

Collaborators: (Lancaster) B. Robinson, S. Jarvis, A. Robson, L. Ponomarenko, N. Tapia-Ruis, S. Rodrigues, G. Adamopoulos.

Colaborators (External) P. Pingue, F. Dinelli (CNR, SNS, Pisa), T. Wang (Sheffield), H. Liu (UCL), P. Smowton, S. Shutts (Cardiff).



Oleg Kolosov's group Bruker European material AFM user meeting 2021



